	Design Specifications	
	OPERATIONAIR: OperationAIR (OperationAIR)	Date: 2020/04/29 14:28:23
	Document: 1	Page: 1 / 12

## 1. TABLE OF CONTENTS

1. Table of Contents.....	1
2. Design specifications .....	1
2.0.1. Features .....	1
2.0.2. Performance .....	4
2.0.3. Safety .....	5
2.0.4. Information provided to the user .....	9

## 2. DESIGN SPECIFICATIONS

### 2.0.1. Features

#### **SPEC-31 Ventilation mode**

Device has a Pressure Controlled Ventilation (PCV) mode

#### **SPEC-46 Usability user-interface**

The user-interface is intuitive, simple, and straightforward to use

#### **SPEC-34 Visualized measurements**

Volume [mL], pressure [cm H<sub>2</sub>O] and flow [L/min] must be visualized on the device's screen in graphs or as numbers

#### **SPEC-150 Instructions For Use**

Instructions For Use (IFU) are supplied with every device and are available for download online

#### **SPEC-42 Compatibility with standard hospital equipment**

The device is compatible with standard hospital equipment according to NEN-EN-ISO 5359:2014/A1:2017 and ISO 5361:2016

#### **SPEC-141 Inspiratory hold**

The device can perform an inspiratory hold to determine Plateau Pressure.  
When the button is pressed, respiration stops at the peak of the inspiration, while measuring the pressure in the system. When the user lets go of the button, respiration instantly starts again

## **SPEC-156 Expiratory Hold**

The device can perform an expiratory hold to determine PEEP.  
When the button is pressed, respiration stops at the end of the expiration, while measuring the pressure in the system. When the user lets go of the button, respiration instantly starts again

## **SPEC-73 Interruption of power supply**

Device is equipped with a UPS to continue ventilation in case of power failure

## **SPEC-60 Breaks**

The device has breaks that shall be used when the device doesn't need to be moved

## **SPEC-27 Flow Sensor**

Device must be equipped with a flow sensor

## **SPEC-153 HEPA filter**

Device can be equipped with a HEPA filter

## **SPEC-25 HME-filter**

Device must be equippable with an ISO 9360-1:2000 certified HME-filter and should be equipped with such a HME filter. These HME filters are disposable and must be changed with every patient.

## **SPEC-166 Ventilator**

The device has a ventilator to mix potential leaking oxygen with surrounding air

### **2.0.1.1. Ventilation Parameters**

## **SPEC-155 Ppeak setting**

The pressure above PEEP can be set by the user in the range 5 to 40 cm H<sub>2</sub>O, in intervals of 1 cm H<sub>2</sub>O

## SPEC-35 PEEP setting

The positive end-expiratory pressure (PEEP) must be adjustable between 5 and 30 cm H<sub>2</sub>O, and the PEEP can be adjusted per 1 cm H<sub>2</sub>O

## SPEC-20 Breathing rate setting

The breathing rate can be set by the user between 5/min and 40/min

## SPEC-38 I:E ratio setting

Inspiratory:expiratory ratio must be adjustable between 2:1 and 1:3

## SPEC-37 FiO<sub>2</sub> setting

The FiO<sub>2</sub> must be adjustable per 5%

### 2.0.1.2. Alarms

## SPEC-6 Alarm for Inspiratory O<sub>2</sub> concentration (FiO<sub>2</sub>)

An alarm for a too high or too low fractional inspiratory oxygen level.

## SPEC-9 Alarm for positive end-expiratory pressure (PEEP)

## SPEC-8 Alarm for Tidal volume (Tv)

Provide an alarm when the tidal volume is outside the accepted range set by the user.

## SPEC-89 Alarm for Peak Pressure (P<sub>peak</sub>)

## SPEC-87 Alarm for Apnea or disconnect

Alarm when inspiratory or expiratory airflow circuit is blocked or disconnected

## SPEC-84 Alarm for empty battery

Device issues a warning when battery is almost empty.

## **SPEC-19 Clear explanation per alarm on GUI**

When an alarm or warning is issued, the device must clearly indicate what alarm is issued

### 2.0.2. Performance

#### **SPEC-95 Error margins: Tidal Volume**

Maximum error between internally and externally measured Tidal Volume is 5%

#### **SPEC-94 Error margins: I:E ratio**

Maximum error between set and measured I:E ratio is 0.05

#### **SPEC-93 Error margins: FiO2**

Maximum error margin between measured and set FiO2 is 5%

#### **SPEC-92 Error margins: breathing rate**

Maximum error between measured and set breathing rates is 1/min

#### **SPEC-90 Error margins: Lung/airway pressure**

Maximum pressure loss between device and patient lungs is 10 cm H2O

#### **SPEC-36 FiO2 range**

Device is able to provide a FiO2 between 21% and 100% at 5% accuracy.

#### **SPEC-26 Tidal Volume**

Device is able to provide tidal volumes between 300 and 700 mL

#### **SPEC-40 Inspiratory Flow**

The inspiratory flow that the device delivers to the patient has a peak value of at least 1.5 L/s

#### **SPEC-118 Reliability**

Under the same settings the device should give the same output.

## SPEC-128 14 day reliability

Device must be able to operate continuously for at least 14 days, preferably more. Specify expected durability in documentation.

### 2.0.3. Safety

## SPEC-29 Peak pressure

Device must keep the peak pressure always below 70 cm H2O

## SPEC-146 Flow sensor shielding

Flow sensor is electrically shielded, all components near the flow sensor are non-conductive or sufficiently grounded

## SPEC-144 Airflow materials

All components through which inspiratory air can flow are

- biocompatible
- oxygen compatible
- clean
- non-porous
- have a high auto-ignition temperature

## SPEC-129 Robustness

The device is able to withstand continuous use and transport

## SPEC-126 Cleaning and disinfection

The device and its materials are able to withstand regular cleaning and disinfection

## SPEC-120 Alarms are easy to understand and intuitive

The alarms are made according to ISO 60601-1-8

## SPEC-113 Impairment of cooling

The medical device must remain safe during the failure of cooling systems, f.e. when ventilation openings are covered.

## **SPEC-103 No air leakage**

The device may not leak air or oxygen

## **SPEC-100 External exhaust outlets**

External exhaust outlets are placed in the bottomplate of the device, where they cannot be blocked by walls or other objects

## **SPEC-99 Temperature of inspired air**

Temperature of the inspired air and the tubes through which the air flows must be below 30 degrees celcius

## **SPEC-85 High Voltage circuit outside casing**

The high voltage circuitry shall be separated from the high O2 network

## **SPEC-83 Testing of the fail-safe valve during production**

The correct functioning of fail-safe valves shall be tested during production

## **SPEC-80 Housing impact test**

The housing must withstand an impact of a solid smooth steel ball, approximately 50 mm in diameter and with a mass of 500 g  $\pm$  25 g, falling freely from 1,3 m height once onto each relevant part of the test sample or swinging like a pendulum, that drops through a vertical distance of 1,3 m, against vertical surfaces. The test is not applied to flat panel displays.

## **SPEC-79 Housing mechanical resistance**

The housing must be able to withstand forces of 250 N without breaking

## **SPEC-72 Leakage current testing**

Since the protective earthing is accessible, leakage current is considered touch current, and the acceptance level for leakage current is <100 uA.

## SPEC-71 Patient leakage current

The values for patient leakage must not exceed 50 uA

## SPEC-70 Touch current

Touch current from or between parts of the medical device within the patient environment must not exceed 100 µA. Leakage current from accessible outer surfaces of the equipment is also considered to be touch current.

## SPEC-67 Acoustic Energy

The acoustic energy emitted by the device shall not exceed 80 dB during worst case normal use

## SPEC-65 Sharp edges

The device shall have no sharp edges that could lead to tearing of protective gloves or the skin

## SPEC-63 Instability unintended force

The device must be able to withstand unintended forces without tipping over.

## SPEC-62 Instability incline surface

Device must not overbalance when placed on an inclined surface at 10 degrees from the horizontal

## SPEC-143 Visibility User Interface

The user interface can be easily used and is clearly visible when wearing protective eyewear

## SPEC-145 Watchdog program

A watchdog program checks every second if every part of the software does what it is supposed to do. When a part of the software malfunctions, the watchdog will instruct a different part to reboot the malfunctioning part.

## SPEC-53 Power Supply

The device shall use a IEC 60601-1 certified mains supply adapter for connection to a 230V AC supply mains.

## **SPEC-32 Expired air**

The expired air must always be separated from inspiratory air

## **SPEC-44 Non-conductive exterior casing**

The exterior casing of the device must be made of a non-conductive material or the material must be grounded

## **SPEC-45 Protective equipment**

The device can be used by users wearing personal protective equipment such as gloves and protective glasses

## **SPEC-23 Protected off switch**

The off switch is protected for accidental pressing

## **SPEC-22 The O2 and air input connectors**

The O2 and air input connectors comply with ISO 18082:2014 and are non-interchangeable

## **SPEC-18 Battery use warning**

When power supply fails and the device starts running on battery power, an alarm is issued

## **SPEC-5 Pressure relief valve**

The device has a mechanical pressure relief valve that opens at 70 cm H2O

## **SPEC-4 Calibrate sensors during assembly**

All pressure sensors and the flow sensor must be calibrated during assembly


## **SPEC-2 Robustness internal components**

All internal components are tightly secured to the bottomplate

## **SPEC-1 Ingress of water**

The exterior casing is spraying water proof (comparable with IP-22)



	Design Specifications	
	OPERATIONAIR: OperationAIR (OperationAIR)	Date: 2020/04/29 14:28:23
	Document: 1	Page: 9 / 12

#### **SPEC-162 Installation Verification**

Installation Verification is performed after assembly:

- Serial numbers of components have to be registered
- The power supply has to be inspected
- The device has to be inspected visually
- Pressure sensors have to be inspected
- The device has to be turned on and off and power down has to be inspected
- Electrical safety has to be inspected
- Sensors, valves and output have to be tested
- Gas system integrity has to be tested.

### 2.0.4. Information provided to the user

#### 2.0.4.1. IFU

#### **SPEC-152 IFU: regulatory requirements**

The medical device must be accompanied by an IFU with referral address of the manufacturer. These documents can be shared electronically or in printing. Additionally the following must be included:

- a list of all items that are part of the medical device,
- instructions for installation, assemblation and modification,
- instructions for cleaning, disinfection, sterilization of each item or the entire device,
- the safety measures that should be applied during installation,
- which parts are suited for the patient environment,
- Measures that must be applied during preventive maintenance,
- the permissible environmental conditions of use of the ME SYSTEM including conditions for transport and storage,
- advice to the RESPONSIBLE ORGANIZATION: – to carry out all adjustment cleaning, sterilization and disinfection PROCEDURES specified therein;
- advice to the RESPONSIBLE ORGANIZATION: – that the assembly of ME SYSTEMS and modifications during the actual service life require evaluation to the requirements of this standard.
- the name of the device;
- the manufacturer of the device;
- the adress of the manufacturer;
- Indications and contra-indications,
- intended use and intended user;
- the clinical benefit for the patient;
- the performance characteristics;
- the degree of accuracy for the parameters of the device;
- how to verify if the right accessories are used;
- undesirable side-effects;

- preperation of the device before use;
- that users must receive a training before using the device;
- information on maintenance and cleaning of the device and accessories;
- information on cleaning and disinfection between the use on different patients;
- all the warnings that will lead to malfunction of the device and/or changes in the working of the device that will lead to risks to the patient. This will include enviromental changes, like diagnostic or therapeutic procedures;
- how the device can safely be disposed;
- publication or printing date with version number;
- that SAE's must be meantioned to the manufacturer;
- how to safely use the software on the device.

## **SPEC-47 IFU: intended use**

IFU contains an explanation of the intended use of the device

## **SPEC-13 IFU: Device not for use during transport**

Warning in IFU that the device may not be used during patient transport

## **SPEC-14 IFU: Check device before use**

Warning in IFU that device functionality must be checked before use

## **SPEC-98 IFU: Do not use when damaged**

Instructions For Use mention that the device must be unplugged when damaged, for example after a fall or when other parts are visibly damaged.

## **SPEC-117 IFU: turn off device before unplugging**

Warning in the IFU that the device has to be turned off on the user interface before cutting the power

## **SPEC-134 IFU: ventilator settings**

The IFU mentions the upper and lower limits of the ventilation settings

## **SPEC-135 IFU: Explanation of the alarms**

The alarms, their causes, and how to deal with them should be clearly explained in the IFU.

## **SPEC-137 IFU: HEPA filter**

The IFU contains instructions about the use of the HEPA-filter

## **SPEC-51 IFU: Instructions for correct and safe maintenance**

The instructions for use shall state instructions for correct and safe maintenance

## **SPEC-138 IFU: alarm boundaries**

The IFU explains how to set the upper and lower alarm boundaries and mentions the alarm limits

## **SPEC-17 Training and clear IFU for the users**

IFU and training materials are developed and validated in collaboration with medical doctors and nurses. IFU warns that the device can only be used by trained users.

## **SPEC-49 IFU: checked by experts**

The Instructions for use shall be reviewed by both internal quality control and external experts

### 2.0.4.2. Label

## **SPEC-151 Label: regulatory requirements**

The label on the device will contain:

- the device name;
- a serial number;
- the manufacturer with address;
- the manufacturing date, as part of the serial number or separately;
- immediate warnings with icons.

All will be written in text or official icons.

## **SPEC-64 Emergency use sticker**

The device has a clear warning sticker with the text that it is a device which is meant for use ONLY in emergencies on COVID-19 patients.

## **SPEC-115 Emergency number on label**

A telephone number is shown on the label that can be called when assistance is required

## **SPEC-133 HEPA filter warning sticker**

Sticker above the expiratory air connector of the device with the text that the HEPA filter must be placed before use

## **SPEC-61 Do not push warning**

The label has a warning to not push (ISO 7010 P017).

## **SPEC-154 Fire warning sticker**

The label contains a warning to keep device away from open fire

## **SPEC-148 Device label**

The label is placed on the device and must be visible from the point of use by the user.

## **SPEC-130 Intuitive for use - training**

Must not require more than 30minutes training for a doctor with some experience of ventilator use.

## **SPEC-48 Vocabulary and Semantics**

Always use the standard nomenclature as defined in BS ISO 19223:2019 regarding lung ventilators and related equipment: Vocabulary and semantics.

## **SPEC-167 Open Source**

All documentation of the AIRone is online available for free. When documentation goes missing or is unreadable, the user can refer to the online database.